


Recombinant Calponin Antibody [clone rCNN1/832] (V3772)

Catalog No.	Formulation	Size
V3772-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	100 ug
V3772-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	20 ug
V3772SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug
V3772IHC-7ML	Prediluted in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide; *For IHC use only*	7 ml


Recombinant **MOUSE MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Rat
Format	Purified
Clonality	Recombinant Mouse Monoclonal
Isotype	Mouse IgG1, kappa
Clone Name	rCNN1/832
Purity	Protein G affinity chromatography
UniProt	P51911
Localization	Cytoplasmic
Applications	Western Blot : 0.5-1ug/ml Flow Cytometry : 1-2ug/10 ⁶ cells Immunohistochemistry (FFPE) : 0.5-1ug/ml for 30 min at RT
Limitations	This recombinant Calponin antibody is available for research use only.



IHC testing of FFPE human uterus stained with recombinant Calponin antibody (clone rCNN1/832). Required HIER: boil tissue sections in 1mM EDTA, pH 7.5-8.5, for 10-20 min followed by cooling at RT for 20 min.



SDS-PAGE analysis of purified, BSA-free recombinant Calponin antibody (clone rCNN1/832) as confirmation of integrity and purity.

Description

Calponin is a calmodulin, F-actin and tropomyosin binding protein which is thought to be involved in the regulation of smooth muscle contraction.

Application Notes

Optimal dilution of the recombinant Calponin antibody should be determined by the researcher.

1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Immunogen

Recombinant human protein was used as the immunogen for the recombinant Calponin antibody.

Storage

Store the recombinant Calponin antibody at 2-8oC (with azide) or aliquot and store at -20oC or colder (without azide).