


Recombinant Creatine kinase B Antibody / CKB / Brain Creatine Kinase [clone rCKBB/8841] (V4373)

Catalog No.	Formulation	Size
V4373-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	100 ug
V4373-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	20 ug
V4373SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug


Recombinant MOUSE MONOCLONAL

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Mouse Monoclonal
Isotype	Mouse IgG1, kappa
Clone Name	rCKBB/8841
Purity	Protein A/G affinity
UniProt	P12277
Localization	Cytoplasm
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml for 30 min at RT Western Blot : 2-4ug/ml
Limitations	This recombinant Creatine kinase B antibody is available for research use only.




Western blot testing of human Y79 cell lysate with recombinant Creatine kinase B antibody. Predicted molecular weight ~43 kDa.



Western blot testing of human HEK293 cell lysate with recombinant Creatine kinase B antibody. Predicted molecular weight ~43 kDa.



IHC staining of FFPE human brain tissue with recombinant Creatine kinase B antibody (clone rCKBB/8841). Inset: PBS used in place of primary Ab (secondary Ab negative control). HIER: boil tissue sections in pH 9 10mM Tris with 1mM EDTA for 20 min and allow to cool before testing.



SDS-PAGE analysis of purified, BSA-free recombinant Creatine kinase B antibody (clone rCKBB/8841) as confirmation of integrity and purity.

Description

Creatine kinases (CK) are a large family of isoenzymes that regulate levels of ATP in subcellular compartments, where they provide ATP at sites of fluctuating energy demand by the transfer of phosphates between creatine and adenine nucleotides. CKs provide the energy of phosphate hydrolysis necessary to drive the normal function of many cellular systems. In cells, the cytosolic CK enzymes consist of two subunits, which can be either B (brain type) or M (muscle type). There are three different isoenzymes: CKMM, CKBB and CKMB. This mAb recognizes the CKBB isoenzyme and does not react with the B subunit in CKMB.

Application Notes

Optimal dilution of the recombinant Creatine kinase B antibody should be determined by the researcher.

Immunogen

Recombinant human full-length protein was used as the immunogen for the recombinant Creatine kinase B antibody.

Storage

Aliquot the recombinant Creatine kinase B antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

