


Recombinant Cytokeratin 17 Antibody [clone RM351] (R20373)

Catalog No.	Formulation	Size
R20373-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM351
Purity	Protein A purified from animal origin-free supernatant
UniProt	Q04695
Localization	Cytoplasmic
Applications	Immunohistochemistry (FFPE) : 1:500-1:1000 Western Blot : 1:1000-1:20000
Limitations	This recombinant Cytokeratin 17 antibody is available for research use only.



Western blot testing of human HeLa cell lysate with recombinant Cytokeratin 17 antibody at 1:10,000 dilution.



IHC staining of FFPE human tonsil tissue with recombinant Cytokeratin 17 antibody at 1:1000.

Description

This antibody reacts to human Cytokeratin 17 (CK-17).

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant Cytokeratin 17 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A peptide corresponding to residues in human Cytokeratin 17 (CK-17) was used as the immunogen for the recombinant Cytokeratin 17 antibody.

Storage

Store the recombinant Cytokeratin 17 antibody at -20oC.