


## Recombinant *Helicobacter pylori* Antibody / *H. pylori* [clone RM400] (R20416)

| Catalog No.  | Formulation                                                      | Size   |
|--------------|------------------------------------------------------------------|--------|
| R20416-0.1ML | Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide | 100 ul |

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

|                    |                                                                                          |
|--------------------|------------------------------------------------------------------------------------------|
| Availability       | 1-3 business days                                                                        |
| Species Reactivity | <i>Helicobacter pylori</i>                                                               |
| Format             | Purified                                                                                 |
| Clonality          | Recombinant Rabbit Monoclonal                                                            |
| Isotype            | Rabbit IgG                                                                               |
| Clone Name         | RM400                                                                                    |
| Purity             | Protein A purified from animal origin-free supernatant                                   |
| UniProt            | N/A                                                                                      |
| Applications       | Immunohistochemistry (FFPE) : 1:100-1:200                                                |
| Limitations        | This recombinant <i>Helicobacter pylori</i> antibody is available for research use only. |



IHC staining of FFPE human stomach tissue with recombinant *Helicobacter pylori* antibody at 1:100.

### Description

This antibody reacts to *Helicobacter pylori*.

### Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant *Helicobacter pylori* antibody may be required due to differences in protocols and secondary/substrate sensitivity.

### Immunogen

Proteins purified from *Helicobacter pylori* was used as the immunogen for the recombinant *Helicobacter pylori* antibody.

## Storage

Store the recombinant *Helicobacter pylori* antibody at -20oC.