


Recombinant Interleukin 2 Antibody [clone IL2/7051R] (V9562)

Catalog No.	Formulation	Size
V9562-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	100 ug
V9562-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	20 ug
V9562SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	IL2/7051R
Purity	Protein A/G affinity
UniProt	P60568
Localization	Secreted
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml
Limitations	This recombinant Interleukin 2 antibody is available for research use only.



IHC staining of FFPE human colon carcinoma tissue with recombinant Interleukin 2 antibody (clone IL2/7051R). HIER: boil tissue sections in pH 9 10mM Tris with 1mM EDTA for 20 min and allow to cool before testing.



SDS-PAGE analysis of purified, BSA-free recombinant Interleukin 2 antibody (clone IL2/7051R) as confirmation of integrity and purity.

Description

Interleukin-2 (IL-2) is a T cell stimulatory cytokine that induces T cell proliferation and NK cell proliferation and activation.

Application Notes

Optimal dilution of the recombinant Interleukin 2 antibody should be determined by the researcher.

Immunogen

A portion of amino acids 21-153 was used as the immunogen for the recombinant Interleukin 2 antibody.

Storage

Aliquot the recombinant Interleukin 2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.