


Recombinant KRT20 Antibody / Cytokeratin 20 [clone KRT20/4380R] (V9717)

Catalog No.	Formulation	Size
V9717-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	100 ug
V9717-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	20 ug
V9717SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG, kappa
Clone Name	KRT20/4380R
Purity	Protein A/G affinity
UniProt	P35900
Localization	Cytoplasm
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml
Limitations	This recombinant KRT20 antibody is available for research use only.



IHC staining of FFPE human colon carcinoma tissue with recombinant KRT20 antibody (clone KRT20/4380R). HIER: boil tissue sections in pH 9 10mM Tris with 1mM EDTA for 20 min and allow to cool before testing.

Description

This MAb recognizes an intermediate filament protein of 46kDa, identified as CK20 (KRT20). KRT20 is abundantly expressed in goblet cells and enterocytes of the gastrointestinal tract. It is a useful marker of pancreatic and colorectal cancer. KRT20 is expressed under normal, hyperplastic and neoplastic conditions. It has been detected in

adenocarcinomas of the colon, stomach and biliary tract. Breast carcinomas are generally non-reactive.

Application Notes

Optimal dilution of the recombinant KRT20 antibody should be determined by the researcher.

Immunogen

A portion of amino acids 324-424 was used as the immunogen for the recombinant KRT20 antibody.

Storage

Aliquot the recombinant KRT20 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.