


## Recombinant Luteinizing Hormone beta Antibody [clone LHb/1612R] (V9347)

| Catalog No.    | Formulation                                                             | Size   |
|----------------|-------------------------------------------------------------------------|--------|
| V9347-100UG    | 0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide | 100 ug |
| V9347-20UG     | 0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide | 20 ug  |
| V9347SAF-100UG | 1 mg/ml in 1X PBS; BSA free, sodium azide free                          | 100 ug |


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

|                    |                                                                                        |
|--------------------|----------------------------------------------------------------------------------------|
| Availability       | 1-3 business days                                                                      |
| Species Reactivity | Human                                                                                  |
| Format             | Purified                                                                               |
| Clonality          | Recombinant Rabbit Monoclonal                                                          |
| Isotype            | Rabbit IgG                                                                             |
| Clone Name         | LHb/1612R                                                                              |
| Purity             | Protein A/G affinity                                                                   |
| UniProt            | P01229                                                                                 |
| Localization       | Cytoplasm, Secreted                                                                    |
| Applications       | Immunohistochemistry (FFPE) : 1-2ug/ml                                                 |
| Limitations        | This recombinant Luteinizing Hormone beta antibody is available for research use only. |



IHC staining of FFPE human pituitary tissue with recombinant Luteinizing Hormone beta antibody (clone LHb/1612R). HIER: boil tissue sections in pH 9 10mM Tris with 1mM EDTA for 20 min and allow to cool before testing.



SDS-PAGE analysis of purified, BSA-free recombinant Luteinizing Hormone beta antibody (clone LHb/1612R) as confirmation of integrity and purity.

## Description

LH has a beta subunit of 121 amino acids (LHB) that confers its specific biologic action and is responsible for interaction with the LH receptor. This beta subunit contains the same amino acids in sequence as the beta subunit of hCG and both stimulate the same receptor; however, the hCG beta subunit contains an additional 24 amino acids and the hormones differ in the composition of their sugar moieties. LH is synthesized and secreted by gonadotrophs in the anterior lobe of the pituitary gland. In concert with the other pituitary gonadotropin follicle-stimulating hormone (FSH), it is necessary for proper reproductive function. In the female, an acute rise of LH levels triggers ovulation. In the male, where LH has also been called Interstitial Cell-Stimulating Hormone (ICSH), it stimulates Leydig cell production of testosterone. LH is a useful marker in classification of pituitary tumors and the study of pituitary disease.

## Application Notes

Optimal dilution of the recombinant Luteinizing Hormone beta antibody should be determined by the researcher.

## Immunogen

Recombinant full-length human LHB protein was used as the immunogen for the recombinant Luteinizing Hormone beta antibody.

## Storage

Aliquot the recombinant Luteinizing Hormone beta antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.