


## Recombinant PR Antibody / Progesterone Receptor [clone RM357] (R20379)

Catalog No.	Formulation	Size
R20379-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM357
Purity	Protein A purified from animal origin-free supernatant
UniProt	P06401
Localization	Nuclear, cytoplasmic
Applications	Immunohistochemistry (FFPE) : 1:100-1:200 Western Blot : 1:1000-1:2000
Limitations	This recombinant PR antibody is available for research use only.



Western blot testing of human T-47D cell lysate with recombinant PR antibody at 1:2000 dilution. Expected molecular weight: 82-94 kDa (PR-A) and 99-120 kDa (PR-B).



IHC staining of FFPE human breast cancer tissue with recombinant PR antibody at 1:100.

## Description

This antibody reacts to human Progesterone Receptor (PR).

## Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant PR antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A peptide corresponding to C-terminus of human Progesterone Receptor was used as the immunogen for the recombinant PR antibody.

## Storage

Store the recombinant PR antibody at -20oC.