


## Recombinant Transferrin Receptor Antibody / TFRC / CD71 [clone RM384] (R20401)

Catalog No.	Formulation	Size
R20401-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM384
Purity	Protein A purified from animal origin-free supernatant
UniProt	P027865
Applications	Immunohistochemistry (FFPE) : 1:250-1:500 Western Blot : 1:1000-1:2000
Limitations	This recombinant Transferrin Receptor antibody is available for research use only.



Western blot testing of human HeLa cell lysate with recombinant Transferrin Receptor antibody at 1:1000 dilution. Predicted molecular weight 85~95 kDa.



IHC staining of FFPE human colon cancer tissue with recombinant Transferrin Receptor antibody at 1:500.



IHC staining of FFPE human placental tissue with recombinant Transferrin Receptor antibody at 1:500.

## Description

This antibody reacts to human TFRC / CD71.

## Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant Transferrin Receptor antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A peptide corresponding to the internal region of the human protein was used as the immunogen for the recombinant Transferrin Receptor antibody.

## Storage

Store the recombinant Transferrin Receptor antibody at -20oC.