


Recombinant TSH Receptor Antibody / B-Chain [clone TSHRB/1607R] (V5501)

Catalog No.	Formulation	Size
V5501-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	100 ug
V5501-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	20 ug
V5501SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG, kappa
Clone Name	TSHRB/1607R
Purity	Protein A/G affinity
UniProt	P16473
Localization	Cytoplasm, Cell membrane
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml
Limitations	This recombinant TSH Receptor antibody is available for research use only.



IHC staining of FFPE human thyroid tissue with recombinant TSH Receptor antibody (clone TSHRB/1607R). HIER: boil tissue sections in pH 9 10mM Tris with 1mM EDTA for 20 min and allow to cool before testing.



SDS-PAGE analysis of purified recombinant TSH Receptor antibody (clone TSHRB/1607R) as confirmation of integrity and purity.

Description

Thyroid-stimulating hormone (TSH, also known as thyrotropin) is a glycoprotein involved in the control of thyroid structure and metabolism, which stimulates the release of the thyroid hormones. TSH is regulated by thyroid hormone (T3) and various retinoid compounds. TSH binds to the thyroid-stimulating hormone receptor (TSHR), which is cleaved into two subunits, A and B, and plays a major role in regulating thyroid function. The third cytoplasmic loop of TSHR has been identified as critical for its role in regulating inositol phosphate and cAMP formation. In Grave s disease, an autoimmune disorder, TSHR is activated by autoantibodies, which may be stimulated by the cleavage of the A and B subunits.

Application Notes

Optimal dilution of the recombinant TSH Receptor antibody should be determined by the researcher.

Immunogen

Recombinant human TSHR B-Chain protein was used as the immunogen for the recombinant TSH Receptor antibody.

Storage

Aliquot the recombinant TSH Receptor antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.