


REST Antibody (R35002)

Catalog No.	Formulation	Size
R35002-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5978
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This REST antibody is available for research use only.



Western blot testing of HEK293 [A] and HeLa [B] lysate with REST antibody at 0.3ug/ml. Predicted molecular weight ~122kDa but routinely observed at ~140kDa.

Description

Additional name(s) for this target protein: RE1-silencing transcription factor

Application Notes

Optimal dilution of the REST antibody should be determined by the researcher.

Immunogen

Amino acids PETPGAEDKGKSSK were used as the immunogen for this REST antibody.

Storage

Aliquot and store the REST antibody at -20oC.