


RICTOR Antibody (F50242)

Catalog No.	Formulation	Size
F50242-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F50242-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	Q6R327
Applications	IHC (Paraffin) : 1:50
Limitations	This RICTOR antibody is available for research use only.



IHC analysis of FFPE human lymph tissue stained with RICTOR antibody

Description

RICTOR and MTOR (FRAP1) are components of a protein complex that integrates nutrient- and growth factor-derived signals to regulate cell growth.

Application Notes

Titration of the RICTOR antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 252-281 from the human protein was used as the immunogen for this RICTOR antibody.

Storage

Aliquot the RICTOR antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.