

RICTOR Antibody (F50243)


Catalog No.	Formulation	Size
F50243-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F50243-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	Q6R327
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:10-1:50
Limitations	This RICTOR antibody is available for research use only.

250
130
95
72

RICTOR antibody western blot analysis in mouse bladder tissue lysate. Predicted molecular weight ~200 kDa



RICTOR antibody immunohistochemistry analysis in formalin fixed and paraffin embedded human kidney tissue.

Description

RICTOR and MTOR (FRAP1) are components of a protein complex that integrates nutrient- and growth factor-derived signals to regulate cell growth.

Application Notes

Titration of the RICTOR antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 1617-1650 from the human protein was used as the immunogen for this RICTOR antibody.

Storage

Aliquot the RICTOR antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.