


RICTOR Antibody (F50245)


Catalog No.	Formulation	Size
F50245-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F50245-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	Q6R327
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:10-1:50
Limitations	This RICTOR antibody is available for research use only.



Western blot analysis of RICTOR antibody and SK-BR-3 lysate.



IHC analysis of FFPE human testis tissue stained with RICTOR antibody

Description

RICTOR is the mTORC2 subunit which regulates cell growth and survival in response to hormonal signals.

Application Notes

Titration of the RICTOR antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 236-267 from the human protein was used as the immunogen for this RICTOR antibody.

Storage

Aliquot the RICTOR antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.