


RIG1 Antibody [ROBO3] (R35750)

Catalog No.	Formulation	Size
R35750-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	64221
Applications	Western Blot : 0.5-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This RIG1 antibody is available for research use only.



Western blot testing of human cerebellum lysate with RIG1 antibody at 0.5ug/ml.
Predicted/observed molecular weight: ~150kDa.

Description

Additional name(s) for this target protein: Retinoblastoma inhibiting gene 1, ROBO3, roundabout 3; ROBO3

Application Notes

Optimal dilution of the RIG1 antibody should be determined by the researcher.

Immunogen

Amino acids RRKQRKELSHYT were used as the immunogen for this RIG1 antibody.

Storage

Aliquot and store the RIG1 antibody at -20oC.