


RIPK4 Antibody (F50790)


Catalog No.	Formulation	Size
F50790-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F50790-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	P57078
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:50-1:100
Limitations	This RIPK4 antibody is available for research use only.



RIPK4 antibody western blot analysis in A549 lysate



IHC analysis of FFPE human breast carcinoma tissue stained with the RIPK4 antibody

Description

The protein encoded by this gene is a serine/threonine protein kinase that interacts with protein kinase C-delta. The encoded protein can also activate NFkappaB and is required for keratinocyte differentiation. This kinase undergoes autophosphorylation.

Application Notes

Titration of the RIPK4 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 802-832 from the human protein was used as the immunogen for this RIPK4 antibody.

Storage

Aliquot the RIPK4 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.