


RNF2 Antibody / RING1B (R35071)

Catalog No.	Formulation	Size
R35071-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6045
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This RNF2 antibody is available for research use only.



Western blot testing of K562 lysate with RNF2 antibody at 0.5 ug/ml. Predicted molecular weight: ~38kDa.

Description

Additional name(s) for this target protein: Ring finger protein 2, RING2, DING, RING1B

Application Notes

Optimal dilution of the RNF2 antibody should be determined by the researcher.

Immunogen

Amino acids PMELYYAPTKEHK were used as the immunogen for this RNF2 antibody.

Storage

Aliquot and store the RNF2 antibody at -20oC.