


RNF39 Antibody (R33213)


Catalog No.	Formulation	Size
R33213-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Rat, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q9H2S5
Gene ID	80352
Localization	Cytoplasmic
Applications	Western Blot : 1-3ug/ml IHC (FFPE) : 10ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This RNF39 antibody is available for research use only.



Western blot testing of human brain lysate with RNF39 antibody at 1ug/ml. Predicted molecular weight: ~46/38 kDa (isoforms 1/2).



IHC testing of FFPE human small intestine tissue with RNF39 antibody at 10ug/ml. Required HIER: steam section in pH6 citrate buffer for 20 min and allow to cool prior to staining.

Description

Additional name(s) for this target protein: Ring finger protein 39

Application Notes

Optimal dilution of the RNF39 antibody should be determined by the researcher.

Immunogen

Amino acids CDPRAPLRIVPAES were used as the immunogen for this RNF39 antibody. This sequence is common to isoforms 1 and 2.

Storage

Aliquot and store the RNF39 antibody at -20oC.