


ROR1 Antibody (R34116)


Catalog No.	Formulation	Size
R34116-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4919
Applications	Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This ROR1 antibody is available for research use only.



IHC testing of FFPE human prostate tissue with ROR1 antibody at 3ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer.



IHC testing of FFPE human prostate tissue with (left) and without (right) ROR1 antibody at 4ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer.

Description

Additional name(s) for this target protein: Receptor tyrosine kinase-like orphan receptor 1

Application Notes

Optimal dilution of the ROR1 antibody should be determined by the researcher.

1. This ROR1 antibody is specific for isoform 1 of the protein.

Immunogen

Amino acids GNKSQKPYKIDSKQA were used as the immunogen for this ROR1 antibody.

Storage

Aliquot and store the ROR1 antibody at -20oC.