


RPE65 Antibody (R34034)

Catalog No.	Formulation	Size
R34034-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Rat
Predicted Reactivity	Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6121
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This RPE65 antibody is available for research use only.



Western blot testing of rat retina lysate with RPE65 antibody at 0.5ug/ml. Predicted molecular weight 60~65kDa.



Western blot testing of rat retina lysate with RPE65 antibody at 0.1ug/ml. Predicted molecular weight 60~65kDa.

Description

Additional name(s) for this target protein: Retinal pigment epithelium-specific protein 65kDa

Application Notes

Optimal dilution of the RPE65 antibody should be determined by the researcher.

Immunogen

Amino acids EVKKNARKAPQPE were used as the immunogen for this RPE65 antibody.

Storage

Aliquot and store the RPE65 antibody at -20oC.