


RPGRIP1L Antibody / Protein fantom (R33497)

Catalog No.	Formulation	Size
R33497-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23322
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This RPGRIP1L antibody is available for research use only.



Western blot testing of human cerebellum lysate with RPGRIP1L antibody at 0.3ug/ml.
Predicted molecular weight: ~151kDa.

Description

Additional name(s) for this target protein: RPGRIP1-like, Fantom, FTM

Application Notes

Optimal dilution of the RPGRIP1L antibody should be determined by the researcher.

Immunogen

Amino acids DKENNKAKRDILK were used as the immunogen for this RPGRIP1L antibody.

Storage

Aliquot and store the RPGRIP1L antibody at -20oC.