


RPS6KA2 Antibody [RSK3] (R34050)


Catalog No.	Formulation	Size
R34050-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Horse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6196
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This RPS6KA2 antibody is available for research use only.



Western blot testing of mouse lung lysate with RPS6KA2 antibody at 0.01ug/ml, observed here at ~75kDa. (1)



Western blot of rat lung lysate with RPS6KA2 antibody at 1ug/ml, observed here at ~75kDa. (1)

Description

Additional name(s) for this target protein: Ribosomal protein S6 kinase, 90kDa, polypeptide 2; RSK3

Application Notes

Optimal dilution of the RPS6KA2 antibody should be determined by the researcher.

1. Molecular weight: predicted at ~83kDa but routinely observed at ~90kDa.

Immunogen

Amino acids EYLSPNQLSRQDVH were used as the immunogen for this RPS6KA2 antibody.

Storage

Aliquot and store the RPS6KA2 antibody at -20oC.