


## RXR alpha Antibody (R33316)

Catalog No.	Formulation	Size
R33316-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6256
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This RXR alpha antibody is available for research use only.



Western blot testing of nuclear HeLa lysate with RXR alpha antibody at 0.3ug/ml.  
Predicted molecular weight: ~51kDa but routinely observed at 54~60kDa.

## Description

Additional name(s) for this target protein: Retinoid X receptor alpha; RXRA

## Application Notes

Optimal dilution of the RXR alpha antibody should be determined by the researcher.

1. This RXR alpha antibody does not react with RXRb or RXRg.

## Immunogen

Amino acids CQVNSSLTSP TGRGSM were used as the immunogen for this RXR alpha antibody.

## Storage

Aliquot and store the RXR alpha antibody at -20oC.

## References (1)