


## RXR beta Antibody (R33317)


Catalog No.	Formulation	Size
R33317-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6257
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This RXR beta antibody is available for research use only.



Western blot testing of A431 lysate and RXR beta antibody at 1ug/ml. Predicted molecular weight: ~57kDa but routinely observed at 57~60kDa.



Western blot of HEK293 lysate overexpressing RXRB tested with RXR beta antibody at 1ug/ml (mock transfection in lane 1). Predicted molecular weight: ~57kDa but routinely observed at 57~60kDa.

## Description

Additional name(s) for this target protein: Retinoid X receptor beta; RXRB

## Application Notes

Optimal dilution of the RXR beta antibody should be determined by the researcher.

1. This RXR beta antibody does not react with RXRa or RXRg.

## Immunogen

Amino acids CRDGMGDSGRDSRSP were used as the immunogen for this RXR beta antibody.

## Storage

Aliquot and store the RXR beta antibody at -20oC.