


S100 Calcium Binding Protein A2 Antibody / S100A2 [clone rS100A2/6482] (V4725)

Catalog No.	Formulation	Size
V4725-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	100 ug
V4725-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	20 ug
V4725SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG2b
Clone Name	rS100A2/6482
Purity	Protein A/G affinity
UniProt	P29034
Localization	Nucleus, Cytoplasm
Applications	ELISA (Order BSA-free Format For Coating) :
Limitations	This S100 Calcium Binding Protein A2 antibody is available for research use only.



SDS-PAGE analysis of purified, BSA-free S100 Calcium Binding Protein A2 antibody (clone rS100A2/6482) as confirmation of integrity and purity.

Description

S100A2 belongs to the S100 family of proteins containing 2 EF-hand calcium-binding motifs. S100A2 may function as a modulator against excess calcium accumulation in normal human mammary epithelial cells and also have a role in suppressing tumor cell growth. This protein may have a tumor suppressor function. Chromosomal rearrangements and

altered expression of this gene have been implicated in breast cancer. Cytoplasmic overexpression may also be of prognostic significance when observed in oral cancer patients, and the S100A2 has also been identified as significantly down regulated in gastric cancer.

Application Notes

Optimal dilution of the S100 Calcium Binding Protein A2 antibody should be determined by the researcher.

Immunogen

Recombinant full-length human protein was used as the immunogen for the S100 Calcium Binding Protein A2 antibody.

Storage

Aliquot the S100 Calcium Binding Protein A2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.