


SALL4 Antibody (F44264)

Catalog No.	Formulation	Size
F44264-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F44264-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	Q9UJQ4
Applications	Western Blot : 1:1000
Limitations	This SALL4 antibody is available for research use only.



SALL4 antibody western blot analysis in 293,HeLa lysate. Isoform A molecular weight: 112~165 kDa, Isoform B molecular weight: 65~95 kDa.

Description

Sall4 is a probable transcription factor. Defects in Sall4 are the cause of Okihiro syndrome; also known as Duane radial ray syndrome (DRRS). It is a disorder characterized by the association of forearm malformations with Duane retraction syndrome. Sall4 is also involved in forelimb and heart development in mice.

Application Notes

Titration of the SALL4 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 1009-1039 from the human protein was used as the immunogen for this SALL4 antibody.

Storage

Aliquot the SALL4 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.