


## SALL4 Antibody (F53118)

Catalog No.	Formulation	Size
F53118-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F53118-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

**Bulk quote request**

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	Q9UJQ4
<b>Applications</b>	Western Blot : 1:1000
<b>Limitations</b>	This SALL4 antibody is available for research use only.



Western blot analysis of lysate from CEM cell line using SALL4 antibody diluted at 1:1000. Isoform A molecular weight: 112~165 kDa, Isoform B molecular weight: 65~95 kDa.

## Description

Transcription factor with a key role in the maintenance and self-renewal of embryonic and hematopoietic stem cells.

## Application Notes

Titration of the SALL4 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

This SALL4 antibody was produced from a rabbit immunized with a KLH conjugated synthetic peptide between 493-526 amino acids from the central region of human SALL4.

## Storage

Aliquot the SALL4 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.