


SAPCD2 Antibody / C9orf140 (F53851)

Catalog No.	Formulation	Size
F53851-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F53851-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q86UD0
Applications	Western Blot : 1:1000-2000
Limitations	This SAPCD2 antibody is available for research use only.



Western blot testing of human 1) HepG2, 2) MCF-7 and 3) PC-3 cell lysate with SAPCD2 antibody at 1:2000. Predicted molecular weight: 43 kDa.

Description

Additional name(s): Suppressor APC domain-containing protein 2; Tumor specificity and mitosis phase-dependent expression protein (TS/MDEP); p42.3

Application Notes

Titration of the SAPCD2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 333-367 from the human protein was used as the immunogen for the SAPCD2 antibody.

Storage

Aliquot the SAPCD2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.