

## SAR1B Antibody (R33465)


| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R33465-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                    |                                                              |
|--------------------|--------------------------------------------------------------|
| Availability       | 1-3 business days                                            |
| Species Reactivity | Human, Mouse, Rat                                            |
| Format             | Antigen affinity purified                                    |
| Clonality          | Polyclonal (goat origin)                                     |
| Isotype            | Goat Ig                                                      |
| Purity             | Antigen affinity                                             |
| Gene ID            | 51128                                                        |
| Applications       | Western Blot : 0.03-0.1ug/ml<br>ELISA (peptide) LOD : 1:8000 |
| Limitations        | This SAR1B antibody is available for research use only.      |



Western blot testing of mouse liver lysate with SAR1B antibody at 0.03ug/ml. Predicted molecular weight: ~22kDa.



Western blot of HEK293 lysate overexpressing SAR1B probed with SAR1B antibody (mock transfection in lane 1). Predicted molecular weight: ~22kDa.

## Description

Additional name(s) for this target protein: SAR1a gene homolog 2

## **Application Notes**

Optimal dilution of the SAR1B antibody should be determined by the researcher.

## **Immunogen**

Amino acids DHERLLESKEELD were used as the immunogen for this SAR1B antibody.

## **Storage**

Aliquot and store the SAR1B antibody at -20oC.