

SARS-CoV-2 Nucleocapsid Antibody (RQ6297)

Catalog No.	Formulation	Size
RQ6297	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	P0DTC9
Applications	ELISA :
Limitations	This SARS-CoV-2 Nucleocapsid antibody is available for research use only.



Description

Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) is an enveloped, positive-sense, single-stranded RNA virus that causes coronavirus disease 2019 (COVID-19). Virus particles include the RNA genetic material and structural proteins needed for invasion of host cells. Once inside the cell the infecting RNA is used to encode structural proteins that make up virus particles, nonstructural proteins that direct virus assembly, transcription, replication and host control and accessory proteins whose function has not been determined. The structural proteins of SARS-CoV-2 include the envelope protein (E), spike or surface glycoprotein (S), membrane protein (M) and the nucleocapsid protein (N). The nucleocapsid phosphoprotein is a structural protein that binds to, protects the viral RNA genome and is involved in packaging the RNA into virus particles. The N protein has been suggested as an antiviral drug target.

Application Notes

Optimal dilution of the SARS-CoV-2 Nucleocapsid antibody should be determined by the researcher.

Immunogen

Amino acids MSDNGPQNQRNAPRITFGGPSDSTGSNQNGERSGARSKQRRPQGLPNNTASWFTALTQHGKEDLKFP
RG QGVPINTNSSPDDQIGYYRRATRRIRGGDGKMKDLSRWY
FY YLGTGPEAGLPYGANKDGIWVATEGALNTPKDHIGTRNPANNAIVLQLPQGTTLPKGFYAEGSRGGSQASSR
SSSR SRNSSRNSTPGSSRGTS ParmagNGGDAALA
LLLLDRLNQLESKMSGKGQQQQGQTVTKKSAAEASKKPRQKRTATKAYNVTQAFGRRGPEQTQGNFGDQELIRQ
GTDY KHWPQIAQFAPSASAFFGMSRIGMEVTPSGTW
LTYTGAIKLDDKDPNFKDQVILLNKHIDAYKTFPTEPKKDKKKKADETQALPQRQKKQQTVTLLPAADLDDFSK
QLQQSM SSADSTQA were used as the immunogen for the SARS-CoV-2 Nucleocapsid antibody.

Storage

After reconstitution, the SARS-CoV-2 Nucleocapsid antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.