


SCARB2 Antibody (R34831)


Catalog No.	Formulation	Size
R34831-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	950
Applications	Western Blot : 0.1-0.3ug/ml Immunohistochemistry (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This SCARB2 antibody is available for research use only.



Western blot testing of human cerebral cortex lysate with SCARB2 antibody at 0.1ug/ml. Predicted molecular weight: ~54/85kDa (unmodified/glycosylated).



IHC staining of FFPE human cerebral cortex with SCARB2 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Scavenger receptor class B, member 2; SR-BII; LIMP2; LIMPII; LGP85; 85 kDa lysosomal sialoglycoprotein scavenger receptor class B, member 2

Application Notes

Optimal dilution of the SCARB2 antibody should be determined by the researcher.

Immunogen

Amino acids NKANIQFGDNGTTIS were used as the immunogen for this SCARB2 antibody.

Storage

Aliquot and store the SCARB2 antibody at -20oC.

References (1)