


SCGN Antibody [clone AFGD-19] (RQ4793)

Catalog No.	Formulation	Size
RQ4793	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human, Mouse, Rat
Format	Purified
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	AFGD-19
Purity	Affinity purified
UniProt	O76038
Applications	Western Blot : 1:500-1:2000
Limitations	This SCGN antibody is available for research use only.



Western blot testing of human fetal brain lysate with SCGN antibody. Predicted molecular weight ~32 kDa.

Description

Secretagogin is expressed at high levels in the pancreatic islets of Langerhans and to a much lesser extent in the gastrointestinal tract (stomach, small intestine and colon), the adrenal medulla and cortex and the thyroid C-cells.

Application Notes

Optimal dilution of the SCGN antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human SCGN was used as the immunogen for the SCGN antibody.

Storage

Store the SCGN antibody at -20oC.