


## SCNN1A Antibody (F54464)


Catalog No.	Formulation	Size
F54464-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F54464-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)


<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Purified
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity purified
<b>UniProt</b>	P37088
<b>Localization</b>	Cytoplasmic, plasma membrane
<b>Applications</b>	Flow Cytometry : 1:25 (1x10e6 cells) Immunofluorescence : 1:25 Immunohistochemistry (FFPE) : 1:25 Western Blot : 1:500-1:2000
<b>Limitations</b>	This SCNN1A antibody is available for research use only.




IHC testing of FFPE human lung carcinoma tissue with SCNN1A antibody. HIER: steam section in pH6 citrate buffer for 20 min and allow to cool prior to staining.




Immunofluorescent staining of human HeLa cells with SCNN1A antibody (green) and DAPI nuclear stain (blue).



Western blot testing of human Daudi cell lysate with SCNN1A antibody. Predicted molecular weight ~76 kDa.



Flow cytometry testing of fixed and permeabilized human WiDr cells with SCNN1A antibody; Blue=isotype control, Green= SCNN1A antibody.



Flow cytometry testing of human WiDr cells with SCNN1A antibody; Blue=isotype control, Green= SCNN1A antibody.

## Description

Nonvoltage-gated, amiloride-sensitive, sodium channels control fluid and electrolyte transport across epithelia in many organs. These channels are heteromeric complexes consisting of 3 subunits: alpha, beta, and gamma. This gene encodes the alpha subunit, and mutations in this gene have been associated with pseudohypoaldosteronism type 1 (PHA1), a rare salt wasting disease resulting from target organ unresponsiveness to mineralocorticoids.

## Application Notes

The stated application concentrations are suggested starting points. Titration of the SCNN1A antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A portion of amino acids 365-391 from the human protein was used as the immunogen for the SCNN1A antibody.

## Storage

Aliquot the SCNN1A antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

