


## Semaphorin 3E Antibody (R34799)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R34799-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                      |                                                                          |
|----------------------|--------------------------------------------------------------------------|
| Availability         | 1-3 business days                                                        |
| Species Reactivity   | Human                                                                    |
| Predicted Reactivity | Dog, Mouse, Pig, Rat                                                     |
| Format               | Antigen affinity purified                                                |
| Clonality            | Polyclonal (goat origin)                                                 |
| Isotype              | Goat Ig                                                                  |
| Purity               | Antigen affinity                                                         |
| Gene ID              | 9723                                                                     |
| Applications         | Immunohistochemistry (FFPE) : 4-6ug/ml<br>ELISA (peptide) LOD : 1:128000 |
| Limitations          | This Semaphorin 3E antibody is available for research use only.          |



IHC testing of FFPE human prostate with Semaphorin 3E antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

### Description

Additional name(s) for this target protein: SEMA3E

### Application Notes

Optimal dilution of the Semaphorin 3E antibody should be determined by the researcher.

### Immunogen

Amino acids NGGRYGTTKDYPDD were used as the immunogen for this Semaphorin 3E antibody.

## Storage

Aliquot and store the Semaphorin 3E antibody at -20oC.