


## Septin 4 Antibody (R35692)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R35692-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                      |                                                            |
|----------------------|------------------------------------------------------------|
| Availability         | 1-3 business days                                          |
| Species Reactivity   | Human, Mouse                                               |
| Predicted Reactivity | Cow, Dog                                                   |
| Format               | Antigen affinity purified                                  |
| Clonality            | Polyclonal (goat origin)                                   |
| Isotype              | Goat Ig                                                    |
| Purity               | Antigen affinity                                           |
| Gene ID              | 5414                                                       |
| Applications         | Western Blot : 1-3ug/ml<br>ELISA (peptide) LOD : 1:64000   |
| Limitations          | This Septin 4 antibody is available for research use only. |



Western blot testing of human heart lysate with Septin 4 antibody at 1ug/ml. Predicted molecular weight: ~53kDa.

## Description

Additional name(s) for this target protein: ARTS (Apoptosis-related protein in the TGF-beta signaling pathway), BRADEION, CE5B3, H5; SEPT4

## Application Notes

Optimal dilution of the Septin 4 antibody should be determined by the researcher.

1. This Septin 4 antibody will detect isoform 1/3/4.

## **Immunogen**

Amino acids RNKLTRESGTD were used as the immunogen for this Septin 4 antibody.

## **Storage**

Aliquot and store the Septin 4 antibody at -20oC.