


## SERPINI1 Antibody / Neuroserpin (F54542)


| Catalog No.   | Formulation                                | Size    |
|---------------|--------------------------------------------|---------|
| F54542-0.4ML  | In 1X PBS, pH 7.4, with 0.09% sodium azide | 0.4 ml  |
| F54542-0.08ML | In 1X PBS, pH 7.4, with 0.09% sodium azide | 0.08 ml |

[Bulk quote request](#)


|                           |                                                                                                                                                   |
|---------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Availability</b>       | 1-3 business days                                                                                                                                 |
| <b>Species Reactivity</b> | Human                                                                                                                                             |
| <b>Format</b>             | Purified                                                                                                                                          |
| <b>Clonality</b>          | Polyclonal (rabbit origin)                                                                                                                        |
| <b>Isotype</b>            | Rabbit Ig                                                                                                                                         |
| <b>Purity</b>             | Antigen affinity purified                                                                                                                         |
| <b>UniProt</b>            | Q99574                                                                                                                                            |
| <b>Applications</b>       | Western Blot : 1:500-1:2000<br>Flow Cytometry : 1:25 (1x10 <sup>6</sup> cells)<br>Immunofluorescence : 1:25<br>Immunohistochemistry (FFPE) : 1:25 |
| <b>Limitations</b>        | This SERPINI1 antibody is available for research use only.                                                                                        |




IHC testing of FFPE human brain tissue with SERPINI1 antibody. HIER: steam section in pH6 citrate buffer for 20 min and allow to cool prior to staining.



Immunofluorescent staining of human HepG2 cells with SERPINI1 antibody (green) and DAPI nuclear stain (blue).



Western blot testing of human SH-SY5Y cell lysate with SERPINI1 antibody. Predicted molecular weight ~46 kDa.



Flow cytometry testing of human HepG2 cells with SERPINI1 antibody; Blue=isotype control, Green= SERPINI1 antibody.

## Description

This gene encodes a member of the serpin superfamily of serine proteinase inhibitors. The protein is primarily secreted by axons in the brain, and preferentially reacts with and inhibits tissue-type plasminogen activator. It is thought to play a role in the regulation of axonal growth and the development of synaptic plasticity. Mutations in this gene result in familial encephalopathy with neuroserpin inclusion bodies (FENIB), which is a dominantly inherited form of familial encephalopathy and epilepsy characterized by the accumulation of mutant neuroserpin polymers. Multiple alternatively spliced variants, encoding the same protein, have been identified.

## Application Notes

The stated application concentrations are suggested starting points. Titration of the SERPINI1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A portion of amino acids 19-45 from the human protein was used as the immunogen for the SERPINI1 antibody.

## Storage

Aliquot the SERPINI1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.