


SETD2 Antibody / HIF1 / KMT3A (R33954)

Catalog No.	Formulation	Size
R33954-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	29072
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This SETD2 antibody is available for research use only.



Western blot testing of mouse fibroblast lysates transiently expressing unrelated and two different siRNA with SETD2 antibody at 0.5ug/ml

Description

Additional name(s) for this target protein: SET domain containing 2

Application Notes

Optimal dilution of the SETD2 antibody should be determined by the researcher.

Immunogen

Amino acids ERDPDKQTQNKE were used as the immunogen for this SETD2 antibody.

Storage

Aliquot and store the SETD2 antibody at -20oC.