


SFRP1 Antibody (R34542)

Catalog No.	Formulation	Size
R34542-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6422
Applications	Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This SFRP1 antibody is available for research use only.



IHC testing of FFPE human kidney with SFRP1 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining. Vesicular staining is seen.

Description

Additional name(s) for this target protein: Secreted frizzled-related protein 1

Application Notes

Optimal dilution of the SFRP1 antibody should be determined by the researcher.

Immunogen

Amino acids KNHECPTFQSVFK were used as the immunogen for this SFRP1 antibody.

Storage

Aliquot and store the SFRP1 antibody at -20oC.