


SH2B3 Antibody / LNK (N-Terminal Region) (F54194)

Catalog No.	Formulation	Size
F54194-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q9UQQ2
Applications	Western Blot : 1:1000-1:2000
Limitations	This SH2B3 antibody is available for research use only.



Western blot testing of human U-251 MG cell lysate with SH2B3 antibody at 1:2000.
Predicted molecular weight ~63 kDa.

Description

Links T-cell receptor activation signal to phospholipase C-gamma-1, GRB2 and phosphatidylinositol 3-kinase.

Application Notes

The stated application concentrations are suggested starting points. Titration of the SH2B3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 67-101 from human SH2B3 was used as the immunogen for the SH2B3 antibody.

Storage

Aliquot the SH2B3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.