


## SIAH1 Antibody / Seven in absentia homolog 1 (R36054)


| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R36054-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                      |                                                                                            |
|----------------------|--------------------------------------------------------------------------------------------|
| Availability         | 1-3 business days                                                                          |
| Species Reactivity   | Human, Mouse, Rat                                                                          |
| Predicted Reactivity | Dog, Cow                                                                                   |
| Format               | Antigen affinity purified                                                                  |
| Clonality            | Polyclonal (goat origin)                                                                   |
| Isotype              | Goat Ig                                                                                    |
| Purity               | Antigen affinity                                                                           |
| Gene ID              | 6477                                                                                       |
| Applications         | Western Blot : 1-3ug/ml<br>Immunoprecipitation : suitable<br>ELISA (peptide) LOD : 1:16000 |
| Limitations          | This SIAH1 antibody is available for research use only.                                    |



Western blot testing of human liver lysate with SIAH1 antibody at 1ug/ml. Predicted molecular weight: ~31kDa.



Western blot of A) mouse brain, B) rat brain, C) mouse liver and D) rat liver lysates with SIAH1 antibody at 1ug/ml. Predicted molecular weight: ~31kDa.

## Description

Additional name(s) for this target protein: Seven in absentia homolog 1

## Application Notes

Optimal dilution of the SIAH1 antibody should be determined by the researcher.

## Immunogen

Amino acids SRQTATALPTGTSKC were used as the immunogen for this SIAH1 antibody.

## Storage

Aliquot and store the SIAH1 antibody at -20oC.