


## Sirt3 Antibody (R35542)

Catalog No.	Formulation	Size
R35542-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23410
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Sirt3 antibody is available for research use only.



Western blot testing of mouse heart lysate with Sirt3 antibody at 0.3ug/ml. The expected ~28kDa band (processed form) and the additional ~15kDa band are both blocked by the immunizing peptide.

## Description

Additional name(s) for this target protein: Sirtuin 3

## Application Notes

Optimal dilution of the Sirt3 antibody should be determined by the researcher.

## Immunogen

Amino acids RASGIPASKLVEAH were used as the immunogen for this Sirt3 antibody.

## Storage

Aliquot and store the Sirt3 antibody at -20oC.