


SLC26A6 Antibody (R35574)

Catalog No.	Formulation	Size
R35574-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	65010
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This SLC26A6 antibody is available for research use only.



Western blot testing of human pancreas lysate with SLC26A6 antibody at 0.1ug/ml.
Predicted molecular weight: ~80kDa.

Description

Additional name(s) for this target protein: Solute carrier family 26, member 6

Application Notes

Optimal dilution of the SLC26A6 antibody should be determined by the researcher.

1. This SLC26A6 antibody will detect isoform 1/2/3/4.

Immunogen

Amino acids RDVAEYSEAKEVR were used as the immunogen for this SLC26A6 antibody.

Storage

Aliquot and store the SLC26A6 antibody at -20oC.