


SLC6A12 Antibody (R33267)

Catalog No.	Formulation	Size
R33267-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6539
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This SLC6A12 antibody is available for research use only.



Western blot testing of human olfactory bulb lysate with SLC6A12 antibody at 1ug/ml. The expected ~70kDa band and the additional ~27kDa band are both blocked by the immunizing peptide.

Description

Additional name(s) for this target protein: Solute carrier family 6 (neurotransmitter transporter, betaine/GABA), member 12

Application Notes

Optimal dilution of the SLC6A12 antibody should be determined by the researcher.

Immunogen

Amino acids CLDGSAGRNFGPSPT were used as the immunogen for this SLC6A12 antibody.

Storage

Aliquot and store the SLC6A12 antibody at -20oC.