


SLC9A3R2 Antibody (isoforms b & e) (R33539)

Catalog No.	Formulation	Size
R33539-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9351
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This SLC9A3R2 antibody is available for research use only.



Western blot testing of Daudi lysate with SLC9A3R2 antibody at 1ug/ml. Predicted molecular weight: ~36kDa.

Description

Additional name(s) for this target protein: Solute carrier family 9 (sodium/hydrogen exchanger), isoform 3 regulator 2

Application Notes

Optimal dilution of the SLC9A3R2 antibody should be determined by the researcher.

1. This SLC9A3R2 antibody is specific for isoforms b and e of the protein.

Immunogen

Amino acids DLPGSDKDTEESGLH were used as the immunogen for this SLC9A3R2 antibody. This sequence is common to isoforms b and e.

Storage

Aliquot and store the SLC9A3R2 antibody at -20oC.