


SMARCA5 Antibody / SNF2H [clone AFDB-19] (RQ4752)

Catalog No.	Formulation	Size
RQ4752	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	AFDB-19
Purity	Affinity purified
UniProt	O60264
Applications	Western Blot : 1:500-1:2000
Limitations	This SMARCA5 antibody is available for research use only.



Western blot testing of human MCF7 cell lysate with SMARCA5 antibody. Predicted molecular weight ~122 kDa.

Description

Helicase that possesses intrinsic ATP-dependent nucleosome-remodeling activity. Complexes containing SMARCA5 are capable of forming ordered nucleosome arrays on chromatin; this may require intact histone H4 tails.

Application Notes

Optimal dilution of the SMARCA5 antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human SNF2H / SMARCA5 was used as the immunogen for the SMARCA5 antibody.

Storage

Store the SMARCA5 antibody at -20oC.