


SNAI3 Antibody (F44977)

Catalog No.	Formulation	Size
F44977-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F44977-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q3KNW1
Applications	Western Blot : 1:1000
Limitations	This SNAI3 antibody is available for research use only.



Description

SNAI3 is a member of the SNAIL gene family, named for the *Drosophila* snail gene, which plays roles in mesodermal formation during embryogenesis.

Application Notes

Titration of the SNAI3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 1-29 from the human protein was used as the immunogen for this SNAI3 antibody.

Storage

Aliquot the SNAI3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.