


SNX15 Antibody (R33747)

Catalog No.	Formulation	Size
R33747-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	29907
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This SNX15 antibody is available for research use only.



Western blot testing of 293 lysate with SNX15 antibody at 0.5ug/ml. Predicted molecular weight: ~38/29kDa (isoforms A/B), observed here at ~45kDa.

Description

Additional name(s) for this target protein: Sorting nexin 15

Application Notes

Optimal dilution of the SNX15 antibody should be determined by the researcher.

Immunogen

Amino acids EEILRLHLSQLPP were used as the immunogen for this SNX15 antibody. This sequence is common to isoforms A and B.

Storage

Aliquot and store the SNX15 antibody at -20oC.