


SORT1 Antibody / Sortilin [clone 17S35] (RQ8887)

Catalog No.	Formulation	Size
RQ8887	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

Recombinant RABBIT MONOCLONAL

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	17S35
Purity	Affinity chromatography
UniProt	Q99523
Applications	Western Blot : 1:500-1:2000
Limitations	This SORT1 antibody is available for research use only.



Western blot testing of human 1) K562, 2) HepG2, 3) U-251 and 4) MCF7 cell lysate with SORT1 antibody. Predicted molecular weight ~92 kDa.