


Sorting nexin 30 Antibody / SNX30 (RQ7753)

Catalog No.	Formulation	Size
RQ7753	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q5VWJ9
Applications	Western Blot : 0.5-1ug/ml Direct ELISA : 0.1-0.5ug/ml
Limitations	This Sorting nexin 30 antibody is available for research use only.



Western blot testing of 1) rat brain and 2) mouse brain tissue lysate with Sorting nexin 30 antibody. Predicted molecular weight ~50 kDa.

Description

Predicted to enable phosphatidylinositol binding activity. Predicted to be involved in protein transport.

Application Notes

Optimal dilution of the Sorting nexin 30 antibody should be determined by the researcher.

Immunogen

E. coli-derived recombinant human protein (amino acids D82-K437) was used as the immunogen for the Sorting nexin 30 antibody.

Storage

After reconstitution, the Sorting nexin 30 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.