


SOX9 Antibody [clone 334CT39.1.4] (F53583)

Catalog No.	Formulation	Size
F53583-0.1ML	In ascites with 0.09% sodium azide	0.1 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Pig
Format	Ascites
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgM
Clone Name	334CT39.1.4
Purity	Ascites
UniProt	P48436
Applications	Western Blot : 1:1000-1:16000
Limitations	This SOX9 antibody is available for research use only.



SOX9 antibody western blot analysis in HepG2 lysate. Predicted molecular weight: 56-65 kDa.

Description

The protein encoded by this gene recognizes the sequence CCTTGAG along with other members of the HMG-box class DNA-binding proteins. It acts during chondrocyte differentiation and, with steroidogenic factor 1, regulates transcription of the anti-Muellerian hormone (AMH) gene. Deficiencies lead to the skeletal malformation syndrome campomelic dysplasia, frequently with sex reversal.

Application Notes

Titration of the SOX9 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 41-70 from the human protein was used as the immunogen for this SOX9 antibody.

Storage

Aliquot the SOX9 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.